



MedGyn

Innovative Products for Ob/Gyn Professionals

- ♀ laminaria
- ♀ Lamicel
- ♀ disposable plastic curettes
- ♀ tubing and collection sets
- ♀ speculums
- ♀ tenaculums
- ♀ dilators
- ♀ forceps
- ♀ needle holders
- ♀ scissors
- ♀ biopsy punches

Established 1976

ISO 9002 CERTIFIED



328 N. EISENHOWER LANE
LOMBARD, ILLINOIS 60148 USA
phone: 630.627.4105
toll free: **800.451.9667**
fax: 630.627.0127

e-mail: medgyn@medgyn.com
web: www.medgyn.com

ob/gyn products for professionals

laminaria

LAMINARIA JAPONICA

This natural cervical dilator is produced from a sea grown plant. Once the plant is dried, it becomes hygroscopic and can swell to approximately four times the dehydrated diameter. This gradual dilation provides less trauma for the patient. Laminaria is available in diameters of 2mm through 10mm with lengths of 60mm and 70mm.

Cat. No.

102	2mm	Extra Small or Extra Thin
102B	2mm	Extra Small or Extra Thin, with Button
103	3mm	Small or Thin
103B	3mm	Small or Thin, with Button
104	4mm	Medium
104B	4mm	Medium, with Button
105	5mm	Large or Thick
105B	5mm	Large or Thick, with Button
105L	5mm	Large or Thick, 70mm length
106	6mm	Extra Large
106B	6mm	Extra Large, with Button
106L	6mm	Extra Large, 70mm length
108	8mm	Jumbo
108L	8mm	Jumbo, 70mm length
110	10mm	Extra Jumbo

GENERAL USE

Laminaria is often utilized when cervical dilation and/or softening of the cervix is desired. Procedures include:

- Treatment of Cervical Stenosis
 - (a) Related to dysmenorrhea
 - (b) Considered a possible cause of infertility
 - (c) Resulting from cauterization or conization
- Induction of labor, drainage of uterine cavity, endometrial biopsy or uterine curette suction by cannula aspiration.
- Placement and removal of intrauterine devices.

MedGyn's Laminaria is clinically proven superior.

- A cervical dilator with more than a century of proven use.
- Provides faster cervical dilation (in only 2-5 hours).
- Expands approximately 4-6 times more.



Comparison of two identical, medium size Laminaria. The one on the right was inserted into cervix at 5:00 pm and removed at 6:30 am the following day (left).

LAMICEL® OSMOTIC CERVICAL DILATOR

Lamicel® is an absorbent polyvinyl acetal sponge, impregnated with less than 500 milligrams of magnesium sulfate (epsom salt). The sponge is compressed to finished size. Once formed, the dilators become highly hygroscopic and capable of increasing their mass up to three times over a period of time, when placed in an environment where they can absorb fluids. Lamicel® OSMOTIC CERVICAL DILATOR, when placed in the endocervical canal, will provide gradual dilation under circumstances found to be tolerable to most patients, and safer and less traumatic than the rapid dilation and attendant risk of employing rigid mechanical dilators.

3mm (average dilation 12mm)
5mm (average dilation 15mm)

cat #113
cat #115



DISPOSABLE VACUUM CURETTES *STERILE* *FLEXIBLE*

The individually packaged sterile, flexible curettes with dual port design, allow the tip to deflect and collapse when pushed against the uterine fundus. This minimizes the risk of perforation. These high quality curettes are available in sizes 4mm through 12mm.

Cat. No.

200	4mm flexible
201	5mm flexible
202	6mm flexible
203	7mm flexible
204	8mm flexible
205	9mm flexible
206	10mm flexible
207	11mm flexible
208	12mm flexible

DISPOSABLE VACUUM CURETTES *STERILE RIGID*

The individually packaged and sterilized, rigid curettes are made from the highest quality clear plastic material to provide strength with a thin wall to provide strength with a thin wall to minimize cervical dilation. These high quality curettes are available straight or curved in sizes 6mm through 16mm.

Cat. No.

Curve

210	6mm curve
211	7mm curve
212	8mm curve
213	9mm curve
214	10mm curve
215	11mm curve
216	12mm curve, large base
217	13mm curve
218	14mm curve
218A	14mm curve w/adaptor
218S	14mm curve, small base
015RC	15mm curve
219	16mm curve

Cat. No.

Straight

220	6mm straight
221	7mm straight
222	8mm straight
223	9mm straight
224	10mm straight
225	11mm straight
226	12mm straight
227	13mm straight
228	14mm straight, large base
228A	14mm straight w/adaptor
228S	14mm straight, small base
015RS	15mm straight
229	16mm straight



Flexible



Rigid Curve



Rigid Straight



Aspiration Kit

ASPIRATION KIT 60cc

The aspiration kit allows accurate, speedy, in-office screening and diagnosis of endometrial cancer and specimen collection. Each kit is individually packaged and sterilized. Available in 4mm, 5mm, and 6mm cannulas.

Cat No.

251	4mm Aspiration Kit
252	5mm Aspiration Kit
253	6mm Aspiration Kit

INTRA-UTERINE SOUNDS

The Intra-Uterine Sound is a single use sound made of white, malleable plastic. Prior to use, the sound can be curved between thumb and finger for following the curve of the cervical canal. With out-patients, the sound can be useful before inserting an IUD or for confirming the position of the uterus. The graduations are in centimeters and allow for the measurement of the uterine cavity.

Cat No.

271



DISPOSABLE SPECULUM

A single use sanitary vaginal speculum. Available in case quantities of 100 pieces.

Cat No.

Case of 100

420	small
425	medium
430	large

PREGNANCY TEST KIT


ONE STEP Hcg Pregnancy Tests

For the Qualitative Detection of Human Chorionic Gonadotropin (Hcg) in urine.

- Rapid result within 5 minutes
- Simple one step procedure
- 25 Min/ML sensitivity for early pregnancy detection
- Easy to read control bands
- Built-in test control
- 18 month shelf life
- Room temperature storage

Cat No.

270



MEDGYN PIPETTE

The MedGyn Suction Pipette has been designed for controlled biopsy of the uterine mucosal lining or menstrual extraction. It has an outer diameter of 2.6mm, which makes it virtually painless to the patient. The MedGyn Pipette is used for cancer screening, endometrial dating, bacterial culturing, menstrual extraction, and more.

Cat No.

262



ENDOSAMPLER w/10cc syringe, 3mm canula / box of 25

Simpler, safer and effective. A unique, sterile, disposable curette for sampling the uterine endometrium.

MedGyn EndoSampler is designed for obtaining clearly differentiated endometrial tissue without anesthesia, in an outpatient setting or as an office procedure. This instrument is designed to improve the efficiency, convenience, and reliability of the endometrial tissue collection and to reduce the physical discomfort of the patient which is usual for the endometrial biopsy. It provides liberal histologic specimens, surgical ease, and is procedurally cost effective for routine use. Except for rare instances, a cervical block, dilator or tenaculum is not required. The MedGyn EndoSampler is recommended for both pre- and post-menopausal women, for diagnostic purposes.


Cat No.

272



EndoSampler

VACUUM CURETTAGE SUPPLIES AND ACCESSORIES



DISPOSABLE COLLECTION SET

Our disposable collection set comes with a 3/8" handle and six feet of tubing individually packaged and sterilized. The swivel handle contains a slip ring for fingertip vacuum control.

Cat No. 231

3/8" Collection Set



DISPOSABLE PLASTIC BOTTLE TOPS

These tops are used only with the 1400cc glass collection bottle. (Cat No. 244) Each top contains a gasket, o-ring and gauze sack, sold ten tops per box.



Cat No. 239 3/8" small conductive
 Cat No. 240 1/2" large conductive
 Cat No. 241 3/8" small non-conductive
 Cat No. 242 1/2" large non-conductive

GLASS COLLECTION BOTTLE

Glass collection bottles with wide-mouth for ease in cleaning.

Cat No. 243 1200cc bottle
 Cat No. 244 1400cc bottle

PERMEABLE GAUZE SACKS

Gauze collection sack readily attaches to the collection bottle for simple and efficient separation of tissue from aspirated fluid. Each package contains ten sacks.

Cat No.

246 (See illustration above)



METAL SWIVEL HANDLES

Reusable handles, with slip ring, for fingertip vacuum control. Size 3/8" holds curettes up to size 14mm with small base. The 1/2" handle holds curettes in sizes 14, 15, and 16mm.

232 3/8" Metal Handle
 235 1/2" Metal Handle
 237 Female fitting (stainless steel)
 238 3/8" Male fitting (stainless steel)
 239 1/2" Male fitting (stainless steel)
 (not shown)

GRAVES

Economy

Cat. No.

- 40-100P (Small, Blades 3/4" x 3")
- 40-110P (Medium, Blades 1 3/8" x 4")
- 40-120P (Large, Blades 1 3/8" x 4 3/4")

GRAVES

Made in Germany

Cat. No.

- 40-090G (Extra-Small Blades 1/2" x 3")
- 40-100G (Small, Blades 3/4" x 3")
- 40-110G (Medium, Blades 1 3/8" x 4")
- 40-120G (Large, Blades 1 3/8" x 4 3/4")
- 40-130G (Extra-Large, Blades 1 1/2" x 6 1/2")



**Graves
Original Model**



**Graves
Improved**

GRAVES IMPROVED

Made in Germany

Cat. No.

- 40-135G (1 1/4" x 4" 45° angle)

GRAVES OPEN SIDE

Made in Germany

Cat. No.

- 40-140G (1 3/8" x 4" medium wide-open angle blades)



**Graves
Open Side**

PEDERSON

Economy

Cat. No.

- 40-200P (Small, Blades 1/2" x 3")
- 40-210P (Medium, Blades 7/8" x 4")
- 40-220P (Large, Blades 1" x 4-3/4")

PEDERSON

Made in Germany

Cat. No.

- 40-180G (Extra Narrow 5/8" x 4")
- 40-190G (Extra Small 5/8" x 2 3/4")
- 40-200G (Small, Blades 1/2" x 3")
- 40-210G (Medium, Blades 7/8" x 4")
- 40-220G (Large, Blades 1" x 4-3/4")
- 40-229P (Extra Narrow, Left Open)
- 40-230P (Extra Narrow, Right Open)



**Pederson
Narrow Blades**

WEISMAN-GRAVES

One side open speculum.

Cat. No.

- 40-150 (1 1/4" x 4, right side open)
- 40-151 (1 1/4" x 4, left side open)



**Weisman-Graves
One Side Open**

KLOPFER

Comfort for the patients, wider opening permits insertion of instruments yet with plenty of space for viewing. Short blade extremely useful for retroverted or antiflex uterus.

Cat. No. 40-326

- (1 3/8" x 3", short blade, wide open)



Klopfer



Short Blade

MOORE-GRAVES

A removal of one inch of stainless steel from the center of the anterior and posterior blades creates this unique speculum. The shorter blades will enable the cervix to be drawn down closer to the introitus, providing better visibility, and a greatly improved working area.

Cat. No.

- 40-310 (short blades 1 x 3 1/4")




Moore-Graves Short Blades

VU-MORE™

Cat. No.

- 40-336 (Medium, blade 1" x 4 1/2",
4cm opening)
40-337 (Large, blade 1" x 5",
4cm opening)



Vu-More™

Graves

VU-MAX™

Cat. No.

- 40-338 (Medium Graves, 1 1/4" x 4",
7cm opening)



Vu-Max™



Graves & Pederson

Vu-More™

Vu-Max™

COLLIN

Cat. No.

- 40-339 (small, 1 1/4" x 3 3/4")
40-340 (medium, 1 3/8" x 4")
40-341 (large, 1 5/8" x 4 1/2")



Collin

CUSCO

Cat. No.

- 40-400 (small, 1 1/4" x 3 1/2"
center screw)
40-401 (medium, 1 3/8" x 4"
center screw)
40-402 (large, 1 3/8" x 4 1/2"
center screw)



Cusco - Center Screw

CUSCO

Cat. No.

- 40-403 (small, 1 1/4" x 3 1/2"
side screw)
40-404 (medium, 1 3/8" x 4"
side screw)
40-405 (large, 1 3/8" x 4 1/2"
side screw)



Cusco - Side Screw

GRUNDMANN

Cat. No.

- 40-350 (small, 1" x 4", wide opened)



Grundmann

AUVARD WEIGHTED

Cat. No.

- 40-320-2 (2 lbs, 1 1/2" x 3 1/4",
slightly angled 75°)
40-320-25 (2 1/2 lbs, 1 1/2" x 3 1/4",
slightly angled 75°)
40-320-45 (2 1/2 lbs, 1 3/8" x 4",
acute angle 45°)
40-320-3 (3 lbs, 1 3/8" x 4",
slightly angled 75°)



Auvard Weighted

STEINER-AUVARD

Cat. No.

- 40-325 (2 1/2 lbs, 1 x 5 1/2",
angled 90°)



Steiner-Auvard


KOGAN ENDOCERVICAL SPECULUM

Cat. No.

- 30-613-1 2-4mm tips, gauge/set screw
30-613-2 4-8mm tips, gauge/set screw



Set Screw




KOGAN ENDOCERVICAL SPECULUM

Cat. No.

- 30-613-3 2-4mm tips, gauge/ratchet
30-613-4 4-8mm tips, gauge/ratchet
30-613-5 8-12mm tips, gauge/ratchet



Ratchet



tenaculums




Schroeder-Braun

SCHROEDER-BRAUN TENACULUM **UTERINE TENACULUM**

9-3/4" Straight (Single Tooth)

Cat. No.
30-700



*Cervical Tenaculum with
Atraumatic Jaws*

CERVICAL TENACULUM **ATRAUMATIC JAWS**

A double row of atraumatic teeth, separated by a groove, gently holds the cervical tissue without crushing or tearing.

Cat. No.
30-610 (1/4" jaw)
30-611 (1/8" jaw)
30-609 (angled, 1/4")



Teale

TEALE TENACULUM **VULSELLUM FORCEPS**

10" Curved, 3 x 4 Teeth, 5mm Wide

Cat. No.
30-712



Skene

SKENE TENACULUM **VULSELLUM FORCEPS**

9-1/2", Curved on Side (Single Tooth)

Cat. No.
30-715

BIERER TENACULUM **CERVICAL TENACULUM**

9", Angled Jaws, 5 x 6 Teeth

Cat. No.
30-612



Bierer

DUPLAY TENACULUM **UTERINE TENACULUM**

11", Double Curved (Single Tooth)

Cat. No.
30-900



Duplay

JACOBS TENACULUM **UTERINE TENACULUM**

8 1/2", Straight,
Serrated Jaws, 2 x 2 Teeth

Cat. No.
30-713

8 1/2", Curved Sideways,
Serrated Jaws, 2 x 2 Teeth

Cat. No.
30-714



Jacobs

SCHROEDER FORCEPS **UTERINE VULSELLUM FORCEPS**

10", Straight (Double Tooth)

Cat. No.
30-727

10", Curved on Flat (Double Tooth)

Cat. No.
30-728

10", Curved Sideways (Double Tooth)

Cat. No.
30-729



*Curved
Sideways*

Curved

Schroeder Straight

JACOBSON TENACULUM **UTERINE TENACULUM**

10", Straight,
Serrated Jaws, 2 x 2 Teeth

Cat. No.
30-952

10", Curved Sideways,
Serrated Jaws, 2 x 2 Teeth

Cat. No.
30-953



Jacobson

ring forceps / needle extenders laminaria inserters

FLETCHER TENACULUM CERVIX HOLDING TENACULUM

9-1/2", Straight


Cat. No.

30-725

9-1/2", Curved

Cat. No.

30-726



Curved

Fletcher Straight

FOERSTER FORCEPS SPONGE FORCEPS

Cat. No. 30-705A (9-1/2", straight, smooth jaw)

Cat. No. 30-705B (9-1/2", curved, smooth jaw)

Cat. No. 30-705 (9-1/2", straight, serrated jaw)

Cat. No. 30-706 (9-1/2", curved, serrated jaw)



Foerster




Laminaria Inserter

LAMINARIA INSERTER

10-1/2", Curved

Cat. No.

30-720




Ramplly with
Serrated Oval Jaws

RAMPLLY TENACULUM UTERINE TENACULUM

10", Straight, Ring Jaw, 12mm

Cat. No.

30-718




Agadi

AGADI FORCEPS 7mm RING FORCEP

9-1/2", Straight, Ring Jaw, 7mm

Cat. No.

30-733



Agadi

AGADI FORCEPS 9mm RING FORCEP

9 3/4", Straight, Ring Jaw, 9mm

Cat. No.

30-732



McMahan

McMAHAN LAMINARIA INSERTER

10-1/2", Curved

Cat. No.

30-721



Concave Tip

NEEDLE ADAPTERS


An affordable means of extending ordinary hypodermic needles when longer needles are required. The paracervical and uterosacral areas can easily be reached with an extended needle.

Cat. No.
32-453 (3")

Cat. No.
32-454 (4")

Cat. No.
32-455 (5")

Cat. No.
32-456 (6")



dilators



Pratt

PRATT UTERINE DILATORS

11-1/2", Double End (set of 8 dilators in Khaki roll, sizes 13-15 to 41-43)

Economy

Cat. No.

30-560P

Made in Germany

Cat. No.


30-560G

PRATT EXTRA LARGE UTERINE DILATORS

11-1/2", Double End (set of 8 dilators in Khaki roll, sizes 45-47 to 77-79)

Cat. No.

30-565



Hegar

HEGAR UTERINE DILATORS

7-1/2", Double End (set of 8 dilators in Khaki roll sizes, 3-4 to 17-18)

Economy


Cat. No.

30-550P

Made in Germany

Cat. No.

30-550G



Hank

HANK UTERINE DILATORS

11", Double End (set of 6 dilators in Khaki roll sizes, 11-12 to 21-22)

Cat. No.

30-500

HERN UTERINE DILATORS

11", Single End
(Set of 21 sizes, 44Fr to 104Fr)

Cat. No.

30-502



Hern

HAWKIN AMBLER UTERINE DILATORS

Single End, Tapered, (Set of 16 dilators sizes, 3/16mm to 18/21mm)

Cat. No.

30-501



*Hawkin
Ambler*

KLEEGMAN CERVICAL DILATORS


The stainless steel Kleegman dilator is patterned after the popular Hank dilators. This instrument is much more narrow and has curved ends and rounded tips of 1.0mm and 2.0mm.

Cat. No.

274



Plastic Denniston



PLASTIC - DENNISTON UTERINE DILATORS

Flexible, Non-disposable, and can be: Boiled, Autoclaved, Gas Sterilized, and Cold Sterilized. (Set of 5 dilators sizes, 5-6mm to 13-14mm)

Cat. No.


30-570

OS CERVICAL DILATOR SET CERVICAL DILATORS

Soft, malleable Teflon material for a more comfortable dilation. Ideal for the post-menopausal patient with stenosis. (Set of 3 dilators, Os Locator, Canal Finder, and the Fundus Dilator)

Cat. No.

273



Fundus Dilator

KLEEGMAN CERVICAL DILATORS

The Teflon (reusable) Kleegman is of a similar size (1.0 and 2.0mm) and function but employs softer and more malleable tips for gentle dilation.

Cat. No.


275

SIMS UTERINE CURETTES

11", Sharp Blades on
Malleable or rigid Shanks

Cat No.

30-209	Sims size #00
30-210	Sims size #0
30-211	Sims size #1
30-212	Sims size #2
30-213	Sims size #3
30-214	Sims size #4
30-215	Sims size #5
30-216	Sims size #6
30-217	Sims size #7
30-218	Sims size #8
30-219	Sims size #9
30-110	Sims size #10
30-112	Sims (set of 12)



Sims



Greene

GREENE UTERINE CURETTES

11", Serrated Blunt Blade on
Malleable Shank

Cat No.

30-012	Size 2
30-014	Size 4
30-016	Size 6



Heaney

HEANEY UTERINE BIOPSY CURETTE

9", Serrated Loop

Cat No.


30-250

THOMAS UTERINE CURETTES

11", Blunt Blades on
Malleable Shanks

Cat No.


30-051	Size 1
30-052	Size 2
30-053	Size 3
30-054	Size 4
30-055	Size 5
30-056	Size 6



Thomas



Bumm



BUMM UTERINE CURETTES

12", Extra Large Sizes on
Malleable Shanks

Cat No.


30-784	Bumm size #1
30-785	Bumm size #2
30-786	Bumm size #3
30-787	Bumm size #4
30-783	set of 4

EVANS UTERINE CURETTES

11", Sharp Blades on Rigid Shanks

Cat No.

30-782	Evans size #0 (4mm)
30-788	Evans size #1 (6mm)
30-789	Evans size #2 (8mm)
30-790	Evans size #3 (10mm)
30-791	Evans size #4 (12mm)
30-792	Evans size #5 (14mm)
30-781	set of 5



Evans

ovum forceps



Bierer

Bierer with Ratchet

BIERER FORCEPS *OVUM FORCEPS*

11", Small, slightly curved, jaws 16mm wide.

Cat. No. 30-760

11", Large, slightly curved for advanced second trimester procedures, jaws 19mm wide.

Cat. No. 30-761

BIERER FORCEPS WITH RACHET *OVUM FORCEPS*

11", Small, slightly curved, jaws 16mm wide.

Cat. No. 30-762

11", Large, slightly curved for advanced second trimester procedures, jaws 19mm wide.

Cat. No. 30-763



*Coarse Serrated
Small (16mm) Oval Jaws*



*Coarse Serrated
Large (19mm) Oval Jaws*



Sopher



*Serrated
Small (12mm) Oval Jaws*



*Serrated
Large (14mm) Oval Jaws*

SOPHER FORCEPS *OVUM FORCEPS*

11", Small, slightly curved, especially useful in second trimester procedures, jaws 12mm wide.

Cat. No.
30-750

11", Large, slightly curved for advanced second trimester procedures, jaws 14mm wide.

Cat. No.
30-755



Bierer

Bierer with Ratchet

BIERER FORCEPS *OVUM FORCEPS*

13", Small, slightly curved, jaws 16mm wide, with gold handles.

Cat. No. 30-764

13", Large, slightly curved for advanced second trimester procedures, jaws 19mm wide, with gold handles.

Cat. No. 30-765

BIERER FORCEPS WITH RACHET *OVUM FORCEPS*

13", Small, slightly curved, jaws 16mm wide, with gold handles.

Cat. No. 30-764R

13", Large, slightly curved for advanced second trimester procedures, jaws 19mm wide, with gold handles.

Cat. No. 30-765R



*Coarse Serrated
Small (16mm) Oval Jaws*



*Coarse Serrated
Large (19mm) Oval Jaws*



Sopher



*Serrated
Large (14mm) Oval Jaws*



*Serrated
Small (12mm) Oval Jaws*



Sopher with Ratchet

SOPHER FORCEPS *OVUM FORCEPS*

13", Small, slightly curved, especially useful in second trimester procedures, jaws 12mm wide, with gold handle.

Cat. No.
30-751

SOPHER FORCEPS *OVUM FORCEPS*

13", Large, slightly curved, for advanced second trimester procedures, jaws 14mm wide, with gold handle.

Cat. No.
30-756

SOPHER FORCEPS *OVUM FORCEPS WITH RATCHET*

13", Small, slightly curved, especially useful in second trimester procedures, jaws 12mm wide, with gold handle.

Cat. No.
30-751R

SOPHER FORCEPS *OVUM FORCEPS WITH RATCHET*

13", Large, slightly curved, for advanced second trimester procedures, jaws 14mm wide, with gold handle.

Cat. No.
30-756R



Sopher "Original"



*Serrated 9mm Oval Jaws
with Concave Tip*



SOPHER "ORIGINAL" FORCEPS *OVUM FORCEPS*

10 1/2", slightly curved, medium teeth, 9mm oval jaw, concave tip to prevent grasping of the fundus,

Cat. No.
30-749



Hern



**Coarse Serrated
Small Oval Jaws**



**Coarse Serrated
Medium Oval Jaws**



**Coarse Serrated
Large Oval Jaws**

HERN FORCEPS EVACUATION FORCEPS

Cat. No.

30-770 (Small 13", Jaw 12mm wide)

Cat. No.

30-771 (Medium 13", Jaw 15mm wide)

Cat. No.

30-772 (Large 13", Jaw 19mm wide)



Hern - Van Lyth



**Serrated
8mm Ring Jaws**



**Serrated
11mm Ring Jaws**

HERN - VAN LYTH FORCEPS EVACUATION FORCEPS

Van Lyth

Cat. No.

30-769 9 1/2", Curved Fine Teeth, 8mm Jaw

Hern-Van Lyth

Cat. No.

30-769A 10", Curved Fine Teeth, 11mm Jaw



Hern - Barrett



**Smooth
16mm Oval Jaws**

HERN - BARRETT FORCEPS PLACENTA FORCEPS

13", Curved, Smooth 16mm Jaw

Cat. No.

30-775



Hern - Patterson



**Coarse Serrated
17mm Ring Jaws**

HERN - PATTERSON FORCEPS EVACUATION FORCEPS

13", Straight with Lock, Coarse
Teeth, 17mm Jaw

Cat. No.

30-768

HERN - KELLY PLACENTA FORCEPS

13", Curved, Medium Teeth,
19mm Jaw

Cat. No.

30-767



**Serrated
19mm Ring Jaws**



Hern-Kelly

ovum forceps



Kelly

KELLY FORCEPS **PLACENTA FORCEPS**

12 1/2", Slightly Curved,
19mm Ring Jaw

Cat. No.
30-766



*Serrated
19mm Ring Jaws*



Boston High Fulcrum

BOSTON FORCEPS **HIGH FULCRUM** **EVACUATION FORCEPS**

13 1/2" Straight, Coarse Teeth,
27mm Jaw

Cat. No.
30-776A



Boston Low Fulcrum

BOSTON FORCEPS **LOW FULCRUM** **EVACUATION FORCEPS**

13 1/2" Straight, Coarse Teeth,
27mm Jaw

Cat. No.
30-776B



*Coarse Serrated
27mm Ring Jaw*



Blumenthal

BLUMENTHAL FORCEPS **EVACUATION FORCEPS**

13", Slightly Curved, Coarse Teeth,
20mm Jaw, concave tip to prevent
grasping of the fundus

Cat. No.
30-777



*Coarse Serrated
20mm Oval
Ring Jaws*



*Concave tip prevents
grasping of the fundus*



Jacobsen

JACOBSEN FORCEPS **EVACUATION FORCEPS**

13", Slightly Curved, Coarse Teeth,
20mm Jaw, Enlarged Ring Handles,
Concave tip to prevent grasping of
the fundus

Cat. No.
30-779



*Coarse Serrated
20mm Oval
Ring Jaws*



*Concave tip prevents
grasping of the fundus*

Enlarged Ring Handles

ovum forceps



*Serrated 13mm
Oval Jaws with
Concave tip to
prevent grasping
of the fundus*

Finks

FINKS FORCEPS **OVUM FORCEPS**

10 1/2", Slightly Curved, Medium
Teeth, 13mm Jaw, concave tip to
prevent grasping of the fundus

Cat. No.
30-736



*Serrated 14mm
Ring Jaw*

McClintock's

McCLINTOCKS FORCEPS **OVUM FORCEPS**

9 1/2", Curved, Fine Teeth,
14mm Ring Jaw

Cat. No.
30-745



13mm x 7.5cm



Cervico Meter

McMAHAN **CERVICO METER**

Angled, used for measuring
the dilation of the cervix.

Cat. No.
30-780



SAENGER **OVUM FORCEPS**

11 1/4" Curved, Oval Jaw with
Inner Diameter Teeth

Cat. No.
60-624

Saenger Curved



SAENGER **OVUM FORCEPS**

10 3/4" Straight with Ratchet,
Oval Jaw with Inner
Diameter Teeth

Cat. No.
60-625

**Saenger Straight
with Ratchet**



SAENGER **OVUM FORCEPS**


10 1/2" Curved with Ratchet,
Oval Jaw with Inner
Diameter Teeth

Cat. No.
60-626

**Saenger Curved
with Ratchet**

OTHER FORCEPS AVAILABLE UPON REQUEST

To Order Call **800.451.9667** MedGyn Products, Inc.



Bozeman Straight

Bozeman Double Curved

BOZEMAN FORCEPS **UTERINE DRESSING FORCEPS**

10-1/2", Curved, stainless steel

Cat. No.

30-710A 10-1/2" straight
30-710B 10-1/2" double curved

ALLIS **TISSUE FORCEPS**

Cat. No.

60-667 4", 4 x 5 Teeth
30-734 5", 4 x 5 Teeth
60-668 6", 4 x 5 Teeth
60-669 6", 5 x 6 Teeth
60-670 7 1/2", 5 x 6 Teeth
60-671 9 1/2", 5 x 6 Teeth
30-735 10", 5 x 6 Teeth
60-673 6", 4 x 5 Teeth, Angular Jaws
60-675 6", 5 x 6 Teeth, Angular Jaws
60-677 7 1/4", 5 x 6 Teeth, Angular Jaws

Cat. No.

Baby Allis
60-680 5 1/2", 4 x 5 Teeth,
Extra Delicate Jaws, 4mm wide
60-681 5", 4 x 5 Teeth, Extra
Delicate Jaws, 2.5mm wide

HALSTED MOSQUITO **HEMOSTATIC FORCEPS**

Straight

Cat. No.

36-156
36-160 (delicate)
36-076 (1 x 2 teeth)


Curved

Cat. No.

36-158
36-162 (delicate)
36-078 (1 x 2 teeth)



Curved 1 x 2 Teeth



Halsted Mosquito



Ochsner

OCSHNER FORCEPS **FORCEPS**

10 1/4", Curved, 1 x 2 Teeth

Cat. No.

30-778

SIMS **UTERINE SOUND, 3mm**

12 1/2", Graduated in Centimeters,
Malleable and Silver Plated.

Cat. No.

30-640A

SIMS **UTERINE SOUND, 2mm**

12 1/2", Graduated in Centimeters,
Malleable and Silver Plated
(2mm dia.)

Cat. No.

30-640B

SIMS **UTERINE SOUND, PROBE**

12 1/2", Plain (**Probe**), Malleable
and Silver Plated

Cat. No.

30-641

SIMPSON **UTERINE SOUND, 3mm**

12 1/2",
Graduated in Inches,

Cat. No.

30-655

12 1/2",
Graduated in Centimeters,

Cat. No.

30-670



Allis



3mm

2mm

Probe

3mm

Sims

Simpson

OTHER FORCEPS AVAILABLE UPON REQUEST



WEBSTER
NEEDLE HOLDER

Cat. No.

60-445 5", Serrated Jaws
60-446 5", Smooth Jaws



MAYO-HEGAR
NEEDLE HOLDER

Cat. No.

60-405-5 5", Serrated Jaws
60-405-6 6", Serrated Jaws
60-405-7 7", Serrated Jaws
60-405-8 8", Serrated Jaws



OLSON-HEGAR
SCISSORS COMBINATION

Cat. No.

60-425 4 3/4", Serrated Jaws
Extra delicate pattern
60-426 4 3/4", Smooth Jaws
Extra delicate pattern
60-427 5 1/2", Serrated Jaws
60-428 6 1/2", Serrated Jaws
60-429 7 1/4", Serrated Jaws



HEANEY
NEEDLE HOLDER

Cat. No.

60-372 8 1/2" Curved Jaws



MATHIEU
NEEDLE HOLDER

Cat. No.

60-380 5 1/2", With Spring Handles
60-390 6 1/4", With Spring Handles
60-391 7 1/2", With Spring Handles



CRILE-WOOD
NEEDLE HOLDER

Cat. No.

60-350-6 6", DeBakey Type
60-350-7 7", DeBakey Type
60-350-8 8", DeBakey Type
60-350-9 9", DeBakey Type
60-350-105 10 1/2", DeBakey Type

OTHER NEEDLE HOLDERS AVAILABLE UPON REQUEST

Scissors



OPERATING SCISSORS STANDARD PATTERN / STRAIGHT

Cat. No.	
25-701	4 1/2", Blunt-Blunt Points
25-703	4 1/2", Sharp-Blunt Points
25-705	4 1/2", Sharp-Sharp Points
25-721	5", Blunt-Blunt Points
25-723	5", Sharp-Blunt Points
25-725	5", Sharp-Sharp Points
25-711	5 1/2", Blunt-Blunt Points
25-713	5 1/2", Sharp-Blunt Points
25-715	5 1/2", Sharp-Sharp Points
25-740	6", Blunt-Blunt Points
25-741	6", Sharp-Blunt Points
25-743	6", Sharp-Sharp Points
25-731	6 1/2", Blunt-Blunt Points
25-733	6 1/2", Sharp-Blunt Points
25-735	6 1/2", Sharp-Sharp Points
25-751	7 1/2", Sharp-Blunt Points



OPERATING SCISSORS STANDARD PATTERN / CURVED

Cat. No.	
25-700	4 1/2", Blunt-Blunt Points
25-702	4 1/2", Sharp-Blunt Points
25-704	4 1/2", Sharp-Sharp Points
25-721	5", Blunt-Blunt Points
25-722	5", Sharp-Blunt Points
25-724	5", Sharp-Sharp Points
25-710	5 1/2", Blunt-Blunt Points
25-712	5 1/2", Sharp-Blunt Points
25-714	5 1/2", Sharp-Sharp Points
25-716	6", Blunt-Blunt Points
25-717	6", Sharp-Blunt Points
25-718	6", Sharp-Sharp Points
25-730	6 1/2", Blunt-Blunt Points
25-732	6 1/2", Sharp-Blunt Points
25-734	6 1/2", Sharp-Sharp Points
25-750	7 1/2", Sharp-Blunt Points



METZENBAUM SCISSORS DISSECTING SCISSORS

Cat. No.	Delicate	Cat. No.	Standard
25-621D	5", Curved	25-621	5", Curved
25-622D	5", Straight	25-622	5", Straight
25-623D	6", Curved	25-623	6", Curved
25-619D	6", Straight	25-619	6", Straight
25-624D	7", Curved	25-624	7", Curved
25-618D	7", Straight	25-618	7", Straight
25-625D	8", Curved	25-625	8", Curved
25-626D	8", Straight	25-626	8", Straight
25-627D	9", Curved	25-627	9", Curved
25-628D	9", Straight	25-628	9", Straight
		25-629	10", Straight
		25-630	10", Curved
		25-631	11", Curved
		25-632	11", Straight
		25-633	12", Straight
		25-634	12", Curved



KELLY-ADSON SCISSORS

Cat. No.	
25-635	6 1/4", Smooth / Straight
25-636	6 1/4", Serrated / Straight
25-637	6 1/4", Smooth / Curved



MAYO SCISSORS DISSECTING SCISSORS

Cat. No.	
25-600	5 1/2", Curved
25-601	5 1/2", Straight
25-602	6 3/4", Curved
25-603	6 3/4", Straight
25-604	9", Curved
25-605	9", Straight
25-606	11", Curved
26-607	11", Straight



DOYEN SCISSORS ABDOMINAL SCISSORS

Cat. No.	
25-000	7", Curved, Smooth
25-005	7", Straight, Smooth
25-010	7", Curved Serrated
25-015	7", Straight Serrated

Scissors



SWEET
45° ANGLE

Cat. No.

25-760 8 1/2", Angled on Flat



LISTER
BANDAGE SCISSORS

Cat. No.

25-150 3 1/2", Extra Fine
25-151 4 1/2", Extra Fine
25-152 5 1/2", Extra Fine
25-153 7 1/4", Extra Fine
25-154 8", With One Large
Finger Ring



LITTAUR
STICH SCISSORS

Cat. No.


25-810 4 1/2", Light Pattern
25-811 5 1/2", Standard Pattern



SPENCER
STICH SCISSORS

Cat. No.


25-850 3 1/2", Delicate



POTTS-SMITH
SCISSORS

Cat. No.


25-671-25 7", Delicate, Angled on
Side 25°, 13mm Blades
25-671-45 7", Delicate, Angled on
Side 45°, 13mm Blades
25-671-60 7", Delicate, Angled on
Side 60°, 13mm Blades
25-670-25 7 1/2", Delicate, Angled on
Side 25°, 20mm Blades
25-670-45 7 1/2", Delicate, Angled on
Side 45°, 20mm Blades
25-670-60 7 1/2", Delicate, Angled on
Side 60°, 20mm Blades



IRIS
SCISSORS

Cat. No.

25-320C 3 1/2", With 20mm Blades,
Curved, Delicate
25-320S 3 1/2", With 20mm Blades,
Straight, Delicate
25-314C 3 3/4", With 14mm Blades,
Curved, Ribbon Type
25-314S 3 3/4", With 14mm Blades,
Straight, Ribbon Type
25-330 4 1/2", Angled on Side
25-335 4 1/2", Curved
25-340 4 1/2", Curved Sideways
25-345 4 1/2", Straight



OTHER SCISSORS AVAILABLE UPON REQUEST

biopsy punches



**TISCHLER-MORGAN™
BIOPSY PUNCH**

Cat. No.

30-110 Bite: 7mm x 3mm



**BABY TISCHLER
BIOPSY PUNCH**

Cat. No.

30-111 Bite: 4.2mm x 2.3mm



**KEVORKIAN
BIOPSY PUNCH**

Cat. No.

30-112 Bite: 8mm x 3mm



**KEVORKIAN-YOUNG
BIOPSY PUNCH**

Cat. No.

30-113 Bite: 8mm x 3mm



**BURKE™
BIOPSY PUNCH**

Cat. No.

30-114 Bite: 6.5mm x 4mm



**EPPENDORFER
BIOPSY PUNCH**

Cat. No.

30-115 Bite: 8mm x 4mm

ROTATING BIOPSY PUNCH

Cat. No.

Shafts

30-113T Tischler Tip
30-113BT Baby Tischler Tip
30-113K Kevorkian Tip
30-113BK Burke Tip
30-113H Handle



TISCHLER

3mm x 7mm Bite



KEVORKIAN

3mm x 8mm Bite




BABY TISCHLER

2.3mm x 4.2mm Bite



BURKE

3mm x 8mm Bite



328 N. EISENHOWER LANE
LOMBARD, ILLINOIS 60148 USA

phone: 630.627.4105

toll free: 800.451.9667

fax: 630.627.0127

e-mail: medgyn@medgyn.com

web: www.medgyn.com