

The new AK 96[®] self-care machine



Confidence in self-care hemodialysis

 **GAMBRO**[®]

Putting patients in control

The compact lightweight AK 96[®] self-care machine allows patients to take an active role in their treatment. The high treatment quality, short training times and cost-efficient operation of the machine make it a solid hemodialysis platform for limited care settings and home care programs.



Easy and safe

With intuitive symbols and clear logical information, the user interface of the AK 96® self-care machine makes handling easy and training less time-consuming.

Multiple settings can be pre-configured or limited, enabling simplified operation and improved safety for patients in self-care dialysis.

The CleanCart® disinfection cartridges, one-hand BPM cuff and battery back-up further enhance the simplicity and safety of the AK 96® self-care machine.

Treatment reliability

The Diascan® function enables a practical and time-saving confirmation of treatment efficacy for both medical professionals and patients. This non-invasive automated real-time clearance monitoring is a valuable complement to blood sampling.

Ultrapure self-care dialysis

The Gambro Hygienic Chain aims at reducing the risk of treatment-related microbiological side effects and inflammatory reactions.

The AK 96® self-care machine equipped with an ultrafilter kit, dry bicarbonate and heat disinfection effectively prevents bacterial growth and biofilm formation. In combination with the single patient Gambro WRO units, this allows the entire fluid circuit to be kept at the same hygienic level.

Integrated self-care HD system

The AK 96® self-care machine and the WRO 300 H unit forms a uniquely integrated self-care HD system allowing all operation from the AK 96® self-care machine. Patient education and training material further complement this offering. The robust design and flexibility of this system make it an ideal partner in any self-care environment.



Convenient

The adjustable remote panel allows the patient to conveniently operate the machine from a sitting or lying position.

Gambro® is a registered trademark of Gambro Lundia AB

AK 96® is a trademark of Gambro Lundia AB registered in the European Community, China and in other countries and filed in the United States

Diascan® is a trademark of Gambro Lundia AB/Gambro Hospal Switzerland Ltd registered in the United States, China and in other countries and filed in the European Community

CleanCart® is a trademark of Gambro Lundia AB/Gambro Hospal Switzerland Ltd, registered in China and in other countries and filed in the United States

Gambro Lundia AB
PO Box 10101
SE-220 10 Lund
Sweden
Phone +46 46 16 90 00
www.gambro.com

